



norden 160r duo superior soft 2.0

Item number: 20082313-00

Detailed information

Item number		Material and surface treatment	
EAN	7392429040839	Color	White
Standards and classifications		Main material	Sanitary acrylic
ETIM	EC011610	Equipment	
EN standard	EN12764, EN61000	Massage system	Superior soft
CE marked	Yes	Heater	Yes
Product dimensions and weight		Outlet drain	Yes
Bathing depth (mm)	430	Drainage pipe	Yes
Height (mm)	670	6 powerjets	Yes
Length (mm)	1570	Back massage	Yes
Width	1150	Foot massage	Yes
Net weight (kg)	115	Mixer	Yes
Water volume (l)	336		Yes
		Packaging dimensions and weight	
		Package height(m)	1,250
		Package width (m)	1,780
		Package depth(m)	0,820
		Package weight(kg)	115,000
		Miscellaneous	
		Model	Offset right

Care instructions

Clean the tub after each bath for long-lasting effect. Use a soft cloth to wipe away the water droplets and prevent limescale deposits. If an acrylic surface gets damaged or scratched, use Clean & Shine Bathroom as this cleans and eliminates minor marks and scratches. For deeper scratches/damage, buff with Restorer. Use our cleaning and disinfection products from SpaCare to keep bacteria away and the pipe system clean.