



motion 140c executive 2.0 + booster

Item number: 20082142-01

Detailed information

Item number		Material and surface treatment	
EAN	7392429041348	Color	White
Standards and classifications		Main material	Sanitary acrylic
ETIM	EC011610	Equipment	
EN standard	EN12764, EN61000	Massage system	Executive 2.0
CE marked	Yes	Compressor	Yes
Product dimensions and weight		Number of air jets	8
Bathing depth (mm)	420	Water cleaning	Yes
Height min. (mm)	675	Numer of water jets	16
Length (mm)	1400	Outlet drain	Yes
Width	1400	Touch panel	Yes
Net weight (kg)	117	Drainage pipe	Yes
Water volume (l)	360	6 powerjets	Yes
Luminaires and electrical components		Back massage	Yes
Number of phases	2	Foot massage	Yes
Incoming voltage (V)	230	Sidemassage	Yes
Power (W)	3150		No
IP rating	IPX5	Packaging dimensions and weight	
		Package height(m)	1,520
		Package width (m)	1,630
		Package depth(m)	0,820
		Package weight(kg)	117,000
		Miscellaneous	
		Model	Corner bathtub

Care instructions

Clean the tub after each bath for long-lasting effect. Use a soft cloth to wipe away the water droplets and prevent limescale deposits. If an acrylic surface gets damaged or scratched, use Clean & Shine Bathroom as this cleans and eliminates minor marks and scratches. For deeper scratches/damage, buff with Restorer. Use our cleaning and disinfection products from SpaCare to keep bacteria away and the pipe system clean.

[Read more on westerbergs.eu](https://www.westerbergs.eu)